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**Preservation of ligament only for export*

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
Advantages and Benefits

- Follow up of 15 years of evolution (2000-2015)
- Excellent wear test performance of different models of global brands of knee joint prostheses market, and demonstrated clinical experience, resulting in increasing the stability and useful life of the implant.
- UHWMPE Polyethylene (Ultra high molecular weight polyethylene), ultra-high strength, biocompatibility, stability, abrasion resistance, and oxidation.
- Cr.Co Alloy: biocompatibility.
- Instrumental Minimally invasive, compact and functional.
- Clinical Studies published by the Department of Orthopedics and Traumatology, School of Medical Sciences of Santa Casa de São Paulo, Pavilion "Fernandinho Simonsen" approved article in Acta Ortop Bras., 2010; 18(6) 310-4.

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
Figura A) Perfil of the femoral implant
 Figura B) Metallic base of the mobile tibial component and central axis for rotation
 Figura C) Polyethylene groove orifice
 Figura D) Free rotating movement, which will be limited by the peripheral structures
 Figura E) Congruence among the components of tibiofemoral joint

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Comparison among wear rates of different knee joint prostheses

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Original Article

TOTAL KNEE ARTHROPLASTY WITH A MOBILE TIBIAL BEARING. MEDIUM-TERM FOLLOW-UP RESULTS

Luiz Gabriel Bezoni Gagliolmetri, Rodrigo de Carmo Costa, Omar Pedro Arês de Camargo, Nilson Roberto Severino, Ricardo de Paula Leite Cary, Victor Marques de Oliveira, Tatiane Alhuza, Roger Araújo

ABSTRACT

Objective: Evaluation of mid-term follow up results of the application of a total complications were one patellar fracture, one distal fracture of the femur, four knee replacements with a mobile tibial bearing design. **Methods:** Ninety four aseptic loosening and four deep joint infections, which required arthroplastic patients (107 knees) were submitted to total knee Arthroplasty, performed with a revision.

Conclusions: With the exception of the cases requiring arthroplastic revision due to the tibial component. The patients were evaluated after a mean follow-up of 32.7 months - standard to septic or aseptic loosening, the authors conclude that the clinical and deviation 21.94 (minimum 24 months and maximum 120 months) through the functional results obtained with Total Knee Replacement with a mobile bearing Knee Society Clinical Rating System (KSSCRS), with a mean outcome of 78.22 component, in a mid-term follow-up, were good.

points. Results: The complications that occurred immediately after or during the surgery included: one wound dehiscence with spontaneous healing, two patellar fractures, one fracture of the medial condyle of the femur, three peroneal nerve outcome. Outcome assessment (health care), impairments, and one sympathetic reflex nervous dystrophy. **Subsequent Postoperative complications.**

Citation: Gagliolmetri LG, Costa RC, Camargo OP, Severino NR, Cary RP, Oliveira VM, et al. Total knee arthroplasty with a mobile tibial bearing. Medium term followup results. Acta Ortop Bras. [online]. 2010;18(6):310-4. Available from URL: <https://www.scielo.br/actb>.